



# Geo5

## K2 CAD / CAM Zirconia

**NEWS**

### Geo5 - K2 CAD / CAM Zirconia

**5 colour gradient Layer**

**5 strength Layer**


**5 translucency Layer**



before



finished glazed work



| layer | height | translucency | strength |
|-------|--------|--------------|----------|
| 1     | 23%    | 47%          | ≥ 800    |
| 2     | 18%    | 46,5%        | ≥ 850    |
| 3     | 18%    | 46%          | ≥ 900    |
| 4     | 18%    | 45%          | ≥ 1000   |
| 5     | 23%    | 44%          | ≥ 1100   |

| Geo5 K2 CAD / CAM Zirconia |          |        |      |      |      |
|----------------------------|----------|--------|------|------|------|
| colour                     | Art-No.  | height |      |      |      |
|                            |          | 14mm   | 16mm | 18mm | 20mm |
| A1                         | 396-90.. | ..14   | ..16 | ..18 | ..20 |
| A2                         | 396-91.. | ..14   | ..16 | ..18 | ..20 |
| A3                         | 396-92.. | ..14   | ..16 | ..18 | ..20 |
| A3,5                       | 396-93.. | ..14   | ..16 | ..18 | ..20 |
| B1                         | 396-94.. | ..14   | ..16 | ..18 | ..20 |
| B2                         | 396-95.. | ..14   | ..16 | ..18 | ..20 |
| C2                         | 396-96.. | ..14   | ..16 | ..18 | ..20 |
| BL1                        | 396-97.. | ..14   | ..16 | ..18 | ..20 |
| BL2                        | 396-98.. | ..14   | ..16 | ..18 | ..20 |